



The Preparation of Low Concentration Hydrogen Sulfide Standards

Gulf Coast Conference

September 12, 2002 Galveston, TX

Paper 050, 12:45 pm Moody F



Robert Benesch, Malik Haouchine, and Tracey Jacksier, Ph.D.

Air Liquide, Chicago Research Center
5230 S. East Ave, Countryside, IL 60525

Previous Focus

- Low level H₂S gas mixtures (<10 ppm) can exhibit reactivity to materials
 - ✓ Reactivity may be observed in flowing systems
 - Long contact times
 - *Low flow rates in tubing*
 - Small contact amounts
 - *Injection of sample into carrier gas stream*
 - ✓ Non-flowing systems
 - Long contact times between gas and material
 - *Storage*

Introduction



6 feet – ¼ inch tubing
Linear velocity 2 cm/sec

Teflon™ - Du Pont

Previous Focus

■ Storage of 100 ppb H₂S in 500 mL sample cylinders

✓ 30 day stability test

- Satisfactory - Sulfinert™, Silcosteel™ -
 - *No experimentally noticeable loss in concentration*
- Unsatisfactory -Aluminum, Stainless Steel, Carbon Steel
 - *Immediate loss of total concentration*

™ Restek Corp.

■ Introducing

Alphatech TM

■ New patent pending, commercial cylinder treatment for specialty gas applications

✓ Provides

- Storage of low level sulfur compounds (≤ 100 ppb) in a commercial compressed gas cylinder
- Stability of H₂S for at least 9 months

Experimental

- To determine the effectiveness of Alphatech™
 - ✓ Cylinders prepared, filled, and analyzed
- Cylinders filled by pressure from source cylinders
 - ✓ 9.9 ppm H₂S
 - ✓ 99.999 % N₂
- Analysis
 - ✓ Utilized 2 different instruments
 - (1) Varian 3800 GC (Initial testing instrument)
 - Sievers 355 SCD (Sulfur Chemiluminescence)
 - Restek 1/16" RT-Sulfur packed column

Experimental

■ Analysis (Con't)

✓ (2) Agilent 6890 GC

- Antek 7090 SCD
- Restek SPB-1 capillary column

✓ Calibration of instruments

- Calibration curves generated within range of testing
 - *Dilution system used to generate four point calibration curve*


✓ Concentration determination

- Averaging four points
- Experimental error (~10 %)

Untreated Cylinder

- Typical compressed gas cylinder – Varied behavior
 - ✓ New aluminum cylinder – Vacuum baking


Untreated Aluminum Cylinder 100 ppb H₂S



Alphatech Stability

- Treated with Alphatech™ process
 - ✓ Three cylinders – 2 filled simultaneously

ALPHATECH™ Stability Test H₂S in Nitrogen



Typical Results

	Initial (calculated)	Final	Avg	RSD
Cylinder 1	79.4	76.4	73.1	5.8%
Cylinder 2	104	111.8	96.6	5.7%
Cylinder 3	104	99.6	99.7	6.4%

→ All cylinders exhibit greater than 9 month stability

■ Adsorption of H₂S onto cylinder walls

✓ Is it occurring in Alphatech™ cylinders?

- Difference between fill concentration and first test concentration
 - *Difference observed is within experimental error*

Fill	Test	% Diff
106.2	111.4	4.90
106.2	105.7	0.47
106.2	108.8	2.45
106.2	111.8	5.27


- Consistency of concentration at various pressures
 - *Increase in concentration at lower pressures if desorption is occurring*

Consistency of Concentration vs. Pressure

■ Cylinder depletion

- ✓ Cylinder filled 97 ppb H_2S in N_2 at production facility using AlphatechTM preparation
- ✓ 3 % difference between 2150 and 63 psia
 - Within experimental error

H_2S in Nitrogen




Operation development

- Currently undergoing a production test with our specialty gas fill plant
 - ✓ Alphatech™ cylinders prepared at specialty gas manufacturing facility
 - ✓ Cylinders filled and analyzed at facility
 - Different size cylinders (16 & 30)
 - Different concentrations (nominal 50 and 100 ppb)

Other Testing


- Application of *Alphatech*TM to other analytes and balance gases

AlphatechTM Stability Test 100 ppb COS in Nitrogen



Other Testing

Alphatech™ Stability Test 50 ppb H₂S in Gaseous Ethylene



Conclusions

- Developed a cylinder treatment / process suitable for storage of low level H₂S compounds
 - ✓ Provide stability for at least 9 months
 - ✓ Currently in operational development phase for production scale up
 - Reliable - Easily repeatable
 - ✓ *Alphatech*TM application to other sulfur analytes and balance gases underway

Production & operations support:

Theo Sadler

Josh Mermelstein

John Kuhn

Ben Duong

■ Booth #130-132

- ✓ Product Specialists
 - Technical Experts