

The Do's and Don'ts in the Analysis of Sulfur for Polyolefin Producers

***Presented as Paper 081 at the
Gulf Coast Conference, Galveston Island, TX
October 22, 2003***

***Benjamin Biela¹, Roger Moore¹,
Robert Benesch², Bruce Talbert² & Tracey Jacksier²***

¹Equistar, Channelview Complex, 8280 Sheldon Rd., Channelview, TX

²Air Liquide, Chicago Research Center, 5230 S. East Ave., Countryside, IL

Agenda

- Introduction
- Ethylene and Propylene
 - ✓ Why is S so important?
 - ✓ What are the stakes?
- Sulfur Analysis
 - ✓ What are the issues?
 - ✓ Analytical Challenges
 - Standard preparation and shelf life
 - How to deliver S to the detector
 - Instrumentation
 - Sample handling
- Summary
- Conclusion

Introduction

- Benjamin Biela, Analytical Specialist
- Equistar Chemicals Analytical Services Laboratory
- Location: Channelview, Texas

Equistar Channelview Analytical Services Lab


- Support 6 olefin plants producing the following products
 - ✓ 11.6 Billion lb/yr Ethylene
 - ✓ 5.0 Billion lb/yr Propylene
- Other products: Butadiene, Benzene, Toluene, Isoprene, Piperylene, Dicyclopentadiene, Alkylate, and Methanol

Sulfur Compounds - Polyolefin Catalyst Poisons

- ppb levels of sulfur compounds present in polyolefin feeds poison catalyst systems and drastically reduce polymer yield
- Reduced yields translate into lost profits (\$\$\$)
 - ✓ Millions of dollars / year
- Typical sulfur impurities
 - ✓ Ethylene
 - Hydrogen Sulfide (H_2S)
 - ✓ Propylene (Polymer Grade)
 - Carbonyl Sulfide (COS)
 - H_2S (Only during startup or unit upset)

Sulfur in Ethylene

- During steady state operation, sulfur in ethylene product is not normally an issue
- A cracked gas compressor surge or trip sometimes results in a caustic wash system upset allowing carbon dioxide and hydrogen sulfide to enter the ethylene splitter and ethylene product



Sulfur in Propylene

- In the cracking furnaces, most (70 >90%) of the inlet COS is converted to H₂S
- Caustic wash towers do not remove a significant amount of COS
- The majority of COS exiting the cracking furnaces will enter the C3 splitter and exit with the propylene product
- A major unit upset or start-up condition could result in a trace amount of H₂S entering the C3 splitter and exiting with the propylene product.

Impact of S on Polymer Yield

Polypropylene Yield Data: Ziegler Natta Catalyst Poison Level for Indicated Effect^{1,2}

Poison	Yield Loss		
	5%	10%	20%
Carbonyl Sulfide	30	60	220
Hydrogen Sulfide	2800	3400	4600

1. *Equistar yield loss data*
2. *Poison level (ppb, wt) for indicated effect*

Specifications

- **Poison levels for Metallocene catalyst: Propylene**
 - ✓ Acceptable maximum sulfur levels
 - Carbon Disulfide (CS₂): 50 ppb
 - Carbonyl Sulfide (COS): 10 ppb
 - Dimethyl Sulfide (C₂H₆S): 1 ppm
- **Manufacturing specifications for polymer grade monomer**
 - ✓ Ethylene
 - Hydrogen Sulfide (H₂S): 50 wt. ppb
 - ✓ Propylene
 - Carbonyl Sulfide (COS): 20 wt. ppb


Analytical Challenges

- If no measurable sulfur response is noted at the analyzer, what does that mean?
 - ✓ The hydrocarbon stream is free from S impurities
 - ✓ “Purified” the hydrocarbon gas stream
 - Reacted out the S

- What are the issues to consider here?

Analytical Challenges

- Preparation of accurate ppb level sulfur standards
- Reasonable shelf life of ppb level sulfur compounds
- Delivery of sulfur compounds to the analytical instrument
 - ✓ Regulator
 - ✓ Tubing
- Analysis
 - ✓ Inlet system
 - ✓ Column



What Exactly is the Challenge?

- Accurate sulfur standards are now available
- Equipment to measure ppb S levels currently exists
- What about the rest of the story?
 - ✓ Regulators
 - ✓ Transfer tubing
 - ✓ GC inlet system
 - ✓ Column

→How do these impact S analysis?

Summary: Storage


- Because of the long contact times associated with storage (minutes to years), any small reactivity of the sulfur compounds with the storage materials becomes critical
 - ✓ A shelf life study was performed with nominal 100 ppb H₂S in N₂ balance mixtures¹
 - Excellent Materials
 - *Silcosteel™, Sulfinert™*
 - » Bonded inert layer to surface of stainless steel
 - *ALPHATECH™ - Air Liquide America*
 - » Cylinder treatment process that yields a shelf life >1 year
 - Poor Materials
 - *Aluminum, Carbon Steel, Stainless Steel*

¹ Benesch, R., Jacksier, T. "The Preparation of Low Concentration Hydrogen Sulfide Standards" GCC 2002, paper 050

Regulators

- Materials not as critical as for storage because of decreased contact times
 - ✓ A regulator comparison study was performed¹ with ppb H₂S in N₂
 - Stainless steel regulators work well once “passivated”
 - *Smaller volume SS regulators are preferred*
 - Materials with brass should be avoided at all costs
 - Sampling flow rates should exceed 100 mL/min

Regulator Study
198 ppb H₂S




→ Signal totally lost with utilization of incorrect regulator

¹ Benesch, R., Haouchine, M., Tabert, B., Jacksier, T. “Low Concentration Hydrogen Sulfide Standards: Is it Possible to Obtain and Regulate at Concentrations below 100 ppb?” Abstract 40 Wednesday Oct 22nd Bluebonnet 9:20 am

Tubing

- Materials are not as critical as for storage because of decreased contact time



Tubing

- Materials – Present in analysis systems
- Sizes 1/4", 1/8", 1/16"

1/4 inch	1/8 inch	1/16 inch
Silcosteel™	Silcosteel™	Silcosteel™
Sulfinert™	Sulfinert™	Sulfinert™
316 SS	316 SS	316 SS
3003 Alum.	3003 Alum.	x
6061 Alum.	x	x
FEP	FEP	x
x	PEEK	PEEK

Procedure

■ Control


- ✓ Prior to testing – control experiments were run without the tubing to determine “control” signals

■ Testing

- ✓ Tubing piece attached and analysis started
 - Pretreatment of tubing
 - *Used “as is” from manufacturer*
 - *Left exposed to air until tested*
 - *No purging after attachment*
 - Analysis finished when signal approached control value

1/4 Inch Tubing

110 ppb Hydrogen Sulfide in Nitrogen
Flow Rate 200 mL/min through 10 feet of 1/4 " Tubing




Surface of 6061 AL was very rough

> 200 minutes before control value was reached

→ Residence time of the gas in the tubing < 1min


1/8 Inch Tubing

110 ppb Hydrogen Sulfide in Nitrogen
Flow Rate 200 mL/min through 10 feet of 1/8 " Tubing



1/16 Inch Tubing


110 ppb Hydrogen Sulfide in Nitrogen
Flow Rate 200 mL/min through 10 feet of 1/16 " Tubing



→ 1/16" SS is not EP and increased time for appearance of S signal maybe related to rougher tubing surface


Step Change in Concentration

- How will a change in the inlet concentration be observed in the outlet concentration?
 - ✓ Will there be a significant time delay in analysis?
 - For 316 SS, up to 200 minute delay per 10 ft of tubing
 - *Lead to process problems*



Step Change In Concentration

Step Change 1/4 " Tubing 350 mL/min : 110 to 150 ppb



- ✓ Silcosteel showed rapid increase
- ✓ 316 SS showed rapid increase once passivated

Passivation of SS


■ Is the passivation of stainless steel permanent?

✓ Three pieces of 316 SS passivated and left to sit under different conditions

1 - Exposed to N₂ for one day

2 - Exposed to air for one day

3 - Exposed to air for one week



✓ After exposure – retested with 100 ppb H₂S in N₂ at 350 mL/min

✓ Long exposures (> 1week) will require repassivation

GC: Valve

■ Loop materials


- ✓ Stainless Steel - requires passivation
- ✓ Silcosteel™, Sulfinert™ - good material

■ Temperature

- ✓ Loops - often heated to prevent adsorption of compounds – how does this effect the reactive sulfurs – 2 concentration levels – 100 ppb and 10,000 ppb
 - Low concentration 100 ppb H₂S in N₂
 - *Silcosteel™ – no noticeable effects from 50 - 225 °C*
 - *Stainless steel – **complete loss of signal** when heated to temperatures above 100 °C, even after passivating*


GC: Sample Loop

- High concentrations 10 ppm
 - *Silcosteel™* - no noticeable effect 50 °C to 225 °C
 - *Stainless Steel* – loss of signal



- On column injection is preferred because of minimal sample contact
- Inlets Tested – materials and temperature effects
 - ✓ Concentration tested – 225 ppb H₂S in N₂
 - Varian 1079
 - *Glass - insert*
 - *Siltek™ insert (Restek)*
 - Agilent Volatiles Interface

Varian 1079 - 225 ppb H₂S in N₂




✓ "Control" determined by using a Silcosteel™ Tee as the inlet


- Volatiles Interface showed no temperature effects (Agilent)

GC: Columns

- Even though the SCD is “sulfur specific”, high concentrations of olefins will produce peaks
- Hydrocarbon peaks must be separated from the S compounds¹
 - ✓ Ethylene (& propylene) and H₂S coelute at temperatures above 0°C
 - ✓ Good separation at elevated temperatures with correct column selection



**100 % dimethyl polysiloxane
split ratio 1:1
-50°C**




¹ Grimberg, A., Benesch, R., Haouchine, M., Richard, S., Coffre, E., Jacksier, T. “Analysis of Sulfur Compounds in Hydrocarbon Matrices”, Pittcon 2002, paper 740-3.

GC: Column

■ Ethylene analysis

- ✓ 100 % Dimethyl polysiloxane
 - often recommended for sulfur analysis
 - Coelution of ethylene and H₂S at non cryogenic conditions
- ✓ Varian Silicaplot
 - Good separation
- ✓ J&W Gas Pro
 - Good separation at temperatures >100°C

Gaspro column
split ratio 4:1
100°C



Summary

- Low concentrations of sulfur will poison Ziegler-Natta and Metallocene catalysts
- Store sulfur compounds in containers with Sulfinert, Silcosteel, or ALPHATECH™ surface treatment technology
- Low wettable surface SS regulators work well
 - ✓ “*Traditional*” SS regulators can be used if passivated and if the flow rate exceeds 100 mL/min
- Tubing should be Sulfinert, Silcosteel, or PEEK
- Column selection is critical for accurate analysis in hydrocarbon matrices
- Sample loop should be coated
- Utilization of “on column” or volatile interface for sample inlet

Conclusion

- Sulfur analysis is not trivial ...the devil is in the details
 - ✓ Correct component selection is the hallmark for a successful analysis
 - Feel sympathy for those who do trace level sulfur analysis
- *Telephone: 281-452-8295*
- *Email: Benjamin.Biela@equistarchem.com*