


Metal-Free LC with Dursan® Bio-Inert Coating

Surface Solutions for Reliable and Repeatable Liquid Chromatography



Overview

A bio-inert flow path is required to manage the complex and reactive sample profiles that are common in today's fast-paced analytical world. Modern LC systems feature an array of exotic materials to achieve bio-inertness, but they are significantly more expensive than instruments constructed with stainless steel flow path components.

Dursan® is a coating for stainless steel that provides equal or better bio-inertness than PEEK along with the robustness of stainless steel. PEEK swells and is not suitable for the pressure conditions required in today's LC. Dursan® is a simple and cost-effective solution for parts requiring bio-inert properties throughout the lab.

Dursan® Specifications

Coating Structure:	Functionalized silica-like coating ($a\text{-SiO}_x\text{:C}_y\text{H}_z$)	
Deposition Process:	Thermal chemical vapor deposition (not plasma-enhanced)	
Temperature:	Deposition	300° to 450°C
	Use	-210°C to 450°C
Substrate:	Compatibility	Stainless Steel, Titanium, Aluminum, Ceramics, More
	Size	Up to 80" (203 cm)
	Geometry	Any shape, including complex geometries
Coating Thickness:	400 - 1600 nm	(Thinner in confined geometries like frits)
Allowable pH Exposure:	0 - 14	
Ideal For:	Frits, Columns, End Fittings, Pump Heads, Valves, Tubing, Needles, and Vials	


Key Features

- Creates a metal-free bio-inert flow path to minimize unwanted analyte interactions and maximize uptime.
- Increases system robustness under extreme salt and pH conditions.
- Inert to nearly all buffers and additives used in HPLC.
- Improves bio-inertness of frits and other difficult components that cannot be treated by other methods.
- Enhances chemical compatibility, even with solvents like tetrahydrofuran (THF) that challenge PEEK.


"The Dursan®-coated columns have so far passed all tests bravely... The results were, as expected, much better than steel columns, but also better than pure PEEK columns."

Performance Data & Benefits

Metal free LC. A powerful metal chelating agent can only be analyzed with a Dursan coated flow path.




Increase corrosion resistance. Eliminate contamination that causes surface activity and instrument downtime.




Improve analytical accuracy. Impurities in this tetracycline sample may have gone unnoticed. Peak tailing is reduced and peak shape improved.

■ Dursan Frit & Column
■ Stainless Steel Frit & Column




Reduce surface fouling from proteins and other biomedica. Increase time between maintenance cycles.



Vaidya, S.V.; Yuan, M.; Narvaez, A.R.; Daghfal, D.; Mattzela, J.; Smith, D. Appl. Surf. Sci. 2016, 364, 896-908

Coat the entire flow path, even frits. A sintered frit was coated, broken in half, and examined via SEM/EDS.



Dark spots are shadowing effects from the rough nature of frit. The entirety of the internal network is coated.



Game-Changing Coatings™
www.SilcoTek.com | +1 (814) 353-1778

Resources

Visit www.SilcoTek.com/learning-center for literature, data, and more.

How to Buy

Visit www.SilcoTek.com/ordering/quote-request for a free custom quote or www.SilcoTek.com/buy-coated-products for stock items.

Patents and Trademarks

SilcoTek® patents and trademarks are the property of SilcoTek Corporation (see <http://www.SilcoTek.com/company-patents-trademarks>). The SilcoTek® registered trademarks used here are registered in the USA and may also be registered in other countries. Information subject to change without notice.

©SilcoTek Corporation. All rights reserved. Printed in the USA.