

Adsorption of ammonia on metal and polymer surfaces


Olavi Vaittinen¹,

Markus Metsälä¹, Stefan Persijn², Markku Vainio^{1,3} and Lauri Halonen¹

¹ University of Helsinki, Finland

² VSL, The Netherlands

³ MIKES, Finland



Variables

Test tube coating: stainless steel 316L, electro-polished SS 316L,
SilcoNert 1000, SilcoNert 2000, Dursan, PFA, FEP, PTFE, PELD, PVDF

NH₃ concentration: 10 ppb – 9 ppm (400 or 9000 ppb)


Flow rate: 0 – 2000 sccm (200 or 1000 sccm)

Temperature: 295 – 333 K (295 K)

Water content: 0 – 1 % (25 ppm)





Gas generation setup





Study protocol

- A) Test tube is flushed with indoor air (≥ 1 h) and pure N_2 (≥ 0.5 h)
- B) Vacuum line and ring-down cavity (except test tube) are exposed to NH_3 (in N_2)
- C) Concentration of ammonia is measured after ~ 1 h
- D) Actual real-time adsorption measurement at 6548.79 cm^{-1} in 3 phases:
 - 1) Ammonia gas flow via by-pass line
 - 2) Ammonia flow switched to go via test tube
 - 3) Slow recovery of ammonia signal



Dursan, 420 ppb





Adsorption on metal / coated surfaces

Metal / coating	Adsorption* (10^{12} molecules/cm ²)	St. deviation (10^{12} molecules/cm ²)
SilcoNert 2000	5.7	0.6
SilcoNert 1000	14.6	0.9
EP SS316L	72	11
Dursan	101	5
SS316L	138	21

*Average of 3 measurements

NH₃ conc = 420 ppb, p (tube) = 176 mbar

Adsorption on polymer surfaces


Polymer	Adsorption* (10^{12} molecules/cm ²)	St. deviation (10^{12} molecules/cm ²)
PVDF	1.0	0.1
PELD	4.4	0.6
PTFE	7.5	1.9
FEP	8.6	0.3
PFA	13.9	1.0

*Average of 3 measurements

NH₃ conc = 8750 ppb, p (tube) = 119 mbar




Langmuir isotherm



Effect of water

SS 316L





Comparison to PTR-MS data*




CRDS:

SS > Dursan >> SN2000 > PFA

PTR-MS:

SS >> Dursan >> SN2000 > PFA

Ammonia transport through different tubings



*SilcoTek Corporation, Sulfinert = SilcoNert 2000

Acknowledgements



The research was carried out with funding by EURAMET and the European Union

The EMRP is jointly funded by the EMRP participating countries within EURAMET and the European Union

More information: **Olavi.Vaittinen@helsinki.fi**

O. Vaittinen *et al*, "Adsorption of gas phase ammonia on treated stainless steel and polymer surfaces,"
submitted to *Measurement Science and Technology*

